

# monsoon

## jet ventilation



### an atmosphere for highest patient comfort



monsoon  
jet ventilation

# monsoon

## jet ventilation



### Technical Specifications

#### Alarm settings:

- PIP Airway Pressure-Alarm: off - 100 mbar
- PP Pause pressure-Alarm: 1 - 100 mbar

#### Monitoring:

- Frequency
- Inspiration time
- FiO<sub>2</sub>
- Airway Pressure; Bypass-Flow
- Alarm limits for: Pause/Airway Pressure, Peak Pressure, Mean Pressure, End Expiratory Pressure, Pause pressure
- TV-/MV monitoring
- etCO<sub>2</sub>

#### Settings:

- Frequency: Jet 1: 4 - 1600/min / Jet 2: 1 - 100/min
- Inspiratory time: 15 - 75 %
- FiO<sub>2</sub>: 21 - 100 %
- Humidification level 0 = 0 %, 8 = 100 %
- Bypass: 0 - 70 l/min
- SIPL: 10 - 40 mbar
- Pressure: 0.1 - 3.5 bar

#### Connections:

- AIR: 4 - 6.5 bar
- O<sub>2</sub>: 4 - 6.5 bar
- Power supply: 100 - 240 VAC, 50/60 Hz
- Power consumption: 100 W

#### Dimensions:

- Size - main unit: 42 x 39 x 19 cm (B x T x H)
- Main unit weight: 14 kg
- Display only: 27 x 21 x 8 cm (B x T x H), Weight: 1.5 k

#### Optional:

- Double Jet Function
- etCO<sub>2</sub>
- Video Camera
- Trolley for accessories



### ACUTRONIC Medical Systems AG

Fabrik im Schiffli | CH-8816 Hirzel

Phone +41 (0)44 729 70 80

Fax +41 (0)44 729 70 81

info@acutronic-medical.ch | www.acutronic-medical.ch